


The determination of Aspirin and caffeine in a mixture by UV spectroscopy



- **Experiment Objective:**
- To appreciate the ability of UV spectrometry in the analysis of a two component mixture.



I. Introduction and theory

Introduction and theory

- UV spectroscopy is generally a non selective technique i.e. presence of more than one substance in the sample would make it difficult to determine the amount of each substance accurately (interference).
- However, in few cases the spectra of two co-present substances in a pharmaceutical formulation might be sufficiently different to allow their accurate determination by UV spectroscopy simultaneously.

I. Introduction and theory

- Aspirin and caffeine represent such an example, where both compounds can be determined simultaneously in their tablet formulation.



I. Introduction and theory

- Aspirin and caffeine are common components of proprietary analgesics. Such products are subjected to strict quality control (QC) procedures to ensure consistency within specified limits.
- The procedures involve various instrumental techniques of which ultraviolet visible spectrometry is of particular importance

I. Introduction and theory

- For a dilute mixture, where the components do not chemically interact in solution and where absorbance curves of each display **clearly defined and separate maxima** quantitative analysis based on Beer- Lambert law is feasible.
- A two component mixture such as aspirin and caffeine in a proprietary analgesic may be analyzed by making absorbance measurements **at two characteristic maxima** (one for each component) and solving for the following pair of simultaneous equations:
 - $A_1 = A_{1a}^{\%} C_a b + A_{1c}^{\%} C_c b$ at $\lambda \max_1$
 - $A_2 = A_{2a}^{\%} C_a b + A_{2c}^{\%} C_c b$ at $\lambda \max_2$

I. Introduction and theory

Note that:

- $A_{\lambda \text{ max1}}$ is a sum of two absorbance
= Absorbance due to first component (aspirin) + Absorbance due to second component (caffeine)
- $A_{\lambda \text{ max1}} =$
($A1\% * B * C$) of first component + ($A1\% * B * c$) of second component

I. Introduction and theory

A1: Mixture
absorbance at
 $\lambda \max_1$

A%1: Caffeine
absorptivity at
 $\lambda \max_1$

$$\bullet A_1 = A_{1a}^{\%} C_a b + A_{1c}^{\%} C_c b \quad \text{at } \lambda \max_1$$

Aspirin
absorptivity at
 $\lambda \max_1$

C_a: Aspirin
concentration

C_c: Caffeine
concentration

$$\bullet A_2 = A_{2a}^{\%} C_a b + A_{2c}^{\%} C_c b \quad \text{at } \lambda \max_2$$

A2: Mixture
absorbance at λ
 \max_2

A%2: Aspirin
absorptivity at
 $\lambda \max_2$

A%2: Caffeine
absorptivity at
 $\lambda \max_2$

II. Practical part



You will need the following Glassware

- (8)25ml Volumetric flasks
- (1) 100ml Volumetric flasks
- (1) 50ml Volumetric flasks
- 0.5 ml volumetric pipette
- 1 ml volumetric pipette

- **Chemicals**

- Solvent:

- methanol

- Analytes:

- Aspirin and caffeine mixed powder
 - Aspirin standard
 - caffeine standard

- **Instrument**

- UV visible spectrophotometer

Procedure

```
graph TD; A[Procedure] --> B[A. Determination of calibration curves]; A --> C[B. Assay of the mixture]
```


A. Determination
of calibration
curves

B. Assay of the
mixture

A. Determination of calibration curves


A. Determination of calibration curves

- **1. Calibration curve for aspirin**
- 1.1 Prepare standard solutions of aspirin




A. Determination of calibration curves

- **2. Calibration curve for caffeine**
- 1.2 Prepare standard solutions of caffeine




UV Scan results

Aspirin spectrum




Caffeine spectrum



A. Determination of calibration curves


- By Now you should have two sets of standard solutions one for aspirin and one for caffeine, and two lambda max one for aspirin and one for caffeine.
- For each set of standard solutions **measure the absorbance at the two lambda max.**



A. Determination of calibration curves

- Using the absorbance values you have you can draw 4 calibration curves, from each you will get an absorptivity value.

Aspirin standard solutions			Caffeine standard solutions		
Standard Number	Absorbance at aspirin lambda max	Absorbance at caffeine lambda max	Standard Number	Absorbance at aspirin lambda max	Absorbance at caffeine lambda max
1 A	0.02	0.12	1 C	0.0024	0.101
2 A	0.31	0.16	2 C	0.0136	0.37
3 A	0.9	0.27	3 C	0.03	0.82
4 A	1.37	0.38	4 C	0.0505	1.25




A. Determination of calibration curves

Calibration curve 1


(Aspirin standard at aspirin lambda max)
Absorptivity (A1% Asp)

Aspirin at 298



Calibration curve 2


(Aspirin standard at caffeine lambda max)
Absorptivity (A2% Asp)



Calibration curve 3

(caffeine standard at aspirin lambda max)
Absorptivity (A1% caff)


Caffeine at 298



Calibration curve 4

(caffeine standard at caffeine lambda max)
Absorptivity = (A2% caff)

caffeine at 272




Attention !
Draw the
calibration
curves to get
the
absorptivity
using the
provided
results by
yourself.
:P

B. Assay of the mixture

B. Assay of the mixture

1. Weigh 315mg of the provided mixture.
2. Dissolve in 50ml methanol and add 20 drops of 4M NaOH.
3. Complete the volume up to 100ml volumetrically.
4. Dilute 1ml up to 50ml with methanol
5. Read the absorbance at the two lambda.



Calculations

- Summary your results:
- Absorptivity values
 - Aspirin At 298= -----
 - Aspirin At 272= -----
 - Aspirin At 298= -----
 - Aspirin At 272= -----
- Mixture absorbance
 - At 298= -----
 - At 272= -----

Calculations

- Substitute your results into the simultaneous equations:

A1: Mixture
absorbance at
 $\lambda \max_1$

A%1: Caffeine
absorptivity at
 $\lambda \max_1$

$$\bullet A_1 = A_{1a}^{\%} C_a b + A_{1c}^{\%} C_c b \quad \text{at } \lambda \max_1$$

Aspirin
absorptivity at
 $\lambda \max_1$

C_a: Aspirin
concentration

C_c: Caffeine
concentration

$$\bullet A_2 = A_{2a}^{\%} C_a b + A_{2c}^{\%} C_c b \quad \text{at } \lambda \max_2$$

A2: Mixture
absorbance at λ
 \max_2

A%2: Aspirin
absorptivity at
 $\lambda \max_2$

A%2: Caffeine
absorptivity at
 $\lambda \max_2$

Calculations

- For example:
- $\lambda \max_1 = 272\text{nm}$ $\lambda \max_2 = 298$
- $A_1 = 0.70$ $A_2 = 0.995$
- $A_{1a}^{\%} = 35.1$ $A_{1c}^{\%} = 595.72$
- $A_{2a}^{\%} = 181.24$ $A_{2c}^{\%} = 24.62$

$$\begin{aligned} 0.70 &= 35.1 C_a + 595.72 C_c & \text{at } \lambda = 272\text{nm} \\ 0.995 &= 181.24 C_a + 24.62 C_c & \text{at } \lambda = 298\text{nm} \end{aligned}$$

You need to solve this pair of equations to get the values of C_c and C_a

Calculations

How to solve simultaneous equation?

To solve by ordinary mathematics follow the link below

- https://www.youtube.com/watch?v=it3vYdV_oyc

www.Math-Mate.com

To solve using excel follow the link below

- <https://www.youtube.com/watch?v=gSNa3fQX0WQ>

Calculations

- Solve the equations to get the values of C_a and C_c
- Multiply C_a and C_c by the dilution factor to get the concentration before dilution.
- Convert the concentration values into amount
 - $\text{Weight (g)} = \text{concentration (g\%)} * \text{volume (ml)}$

Calculations

- The solution for previous example :
- **By solving simultaneous equation using excel:**
 - $C_a = 0.0053\text{g\%}$
 - $C_c = 0.000859\text{g\%}$
- **The concentration before dilution:**
 - $C_a = 0.00537 * 50 = 0.2685 \text{ g\%}$
 - $C_c = 0.000859 * 50 = 0.043 \text{ g\%}$
- **Amount in g**
 - Aspirin = $0.2685 * 100 / 100 = 0.2685\text{g}$ (268.5mg)
 - Caffiene = $0.043\text{g} * 100 / 100 = 0.448 \text{ gm}$ (44.8 mg)



Ca	Cc		constant
35.1	595.27		0.7
181.24	24.62		0.995
-0.00023	0.005562	Ca =	0.005373
0.001693	-0.00033	Cc =	0.000859